

Automated Haematology Analyser
XN-L Series

Leading Haematology
for Better Patient Care



Upgrade your 3PD to 6PD with Sysmex

Why 5PD?

3PD Analyser

Neutrophils

Mixed

Lymphocytes

5PD Analyser

Neutrophils

Eosinophils

Basophils


Monocytes

Lymphocytes

Five part differential analysers count monocytes, eosinophils and basophils separately rather than as a mixed population in three part differential analysers. These individual WBC parameters reveal more valuable information to clinicians to support diagnostic and treatment decisions.


Why Sysmex 6PD?

3PD Analyser



Sysmex 6PD Analyser

added-value
IG¹



Examples of WBC Flags:

WBC Abn Scattergram
Neutropenia
Neutrophilia
Lymphopenia
Lymphocytosis
Monocytosis
Eosinophilia
Leukocytopenia
Leukocytosis
Blasts/Abn Lympho?
Left Shift?
Atypical Lympho?
NRBC?

Benefits of Sysmex 6PD Analyser

1) Assessment is not solely based on cell size

Unlike three part differential and some five part differential analysers, fluorescent flow cytometry measures not only the cell size, but also the intracellular information and nucleic acid content. This produces a highly accurate differential count in EDTA blood samples even as the cell size changes during normal storage.

2) Identification of immature cells

Identification of immature cells is possible with the XN-L as immature cells have a higher nucleic acid content. This has made the generation of six part differential, the IG count possible. The precision of an automated IG count increases laboratory efficiency by reducing manual counts.

3) Superior flagging system for abnormal cells

As compared to a flagging system based entirely on cell size on a three part differential analyser, the XN-L provides more detailed and specific flagging for abnormal cells.

Introducing XN-L Series

XN-L Series is the latest compact fully-automated 6-part differential haematology analyser from Sysmex.

It is designed to meet today's laboratory needs by providing enhanced clinical values that only high-end models were previously able to provide; delivering improved operational efficiency in the laboratories.

XN-L Series is available in four models with different aspiration modes.



XN-550

Sampler that allows both open and closed tube analysis



XN-450

Closed tube analysis



XN-330 / XN-350

Open tube analysis

A completely new compact 6-part differential haematology analyser that caters for what matters most in your lab.

1 Ease of use & peace of mind

Multi-coloured touch screen display embedded with intuitive graphical user interface to support your routine tasks.



A reagent management system that calculates remaining tests and expiry dates to allow more efficient tracking of reagent usage.



2 Space saving

Compact design with an in-built IPU that fits easily to any laboratory bench or table with a smaller footprint.



3 Flexibility

Only **25µL** of aspiration volume

Requires only 25 µl of aspiration volume for whole blood mode and low WBC mode.

Auto-diluent dispensing function

Auto-diluent dispensing function is available for dilution of samples

4 Improved workflow and TAT

Up to **70** samples/hour*

Better TAT with one of the highest throughput compact 6PD analysers in the market.

Sampler available for a truly walk-away system**




Left drawer



Right drawer

Sampler comes with two drawers to allow continuous loading of samples

Integrated auto Repeat, Rerun and Reflex** of sample




*Only available on XN-350 / XN-450 / XN-550

**Only available on XN-550

More than just CBC + DIFF

Core Technology of XN-L Series

XN-L Series utilises the same core technology as high-end haematology analysers - the proven and reliable laser flow cytometry - to count and analyse cells. After the cells are irradiated by the laser beam, the forward scattered light (FSC), side scattered light (SSC) and side fluorescent light (SFL) of the cells are analysed. The three signals are used to differentiate and count cells with the help of unique digital technology and algorithms.



Standard application

added-value
IG^{1,2}

WDF Scattergram

Added value:
A true 6-part differential count (IG):

- Immature granulocytes (IG) with every differential analysis aids in the early prediction of infection and inflammation.
- **No additional reagents required.**

Optional applications*

added-value

Low WBC


Added value:

Analysis of low WBC samples in L-WBC mode:

- Samples with low white blood cell count can be measured accurately by counting twice as many cells; providing reliable results that aid in chemotherapy monitoring.

added-value


Body Fluid^{3,4}



Added value:

Fully-automated body fluid analysis in BF mode:


- 2-part differential body fluid analysis includes MN (mononuclear) and PMN (polymorphonuclear) cell population to aid in the distinction between viral and bacterial infection.



- **No additional reagents required.**
- No special sample preparation required.

added-value

RET⁵




Added value:

Indices of erythropoiesis (RET, RET-He, IRF):

- RET-He (reticulocytes haemoglobin) and IRF (immature reticulocytes fraction) aids in monitoring of RBC production.
- RET-He (reticulocytes haemoglobin) aids in differentiation between functional and classical iron deficiency and monitoring of EPO and/or IV iron therapy.

added-value

IPF^{6,7,8}



Added value:

Indices of thrombopoiesis (IPF):


- IPF (Immature platelet fraction) aids in differential diagnosis of thrombocytopenic disorders and is an early predictor of platelet recovery.

*Only available on XN-350 / XN-450 / XN-550

What's more

Extended Connectivity

- XN-L supports extended connectivity to meet a laboratory's various connectivity requirements.



Standardisation across Sysmex platforms

- Compatibility with the XN-Series through common reagents, controls and technology allows for consistency across laboratory operations.

1 Common reagents and QC materials




Reagents and QC materials utilised with the XN-Series can be used on the XN-L Series.

2 Uniform user interface



The XN-L analysers use the same user interface as the XN-Series.

3 Identical core technology and clinical parameters



Enhanced clinical parameters makes the XN-L Series the perfect backup for 24/7 laboratories.

Specifications



XN-330



XN-350



XN-450



XN-550

Principles & Technologies

Fluorescent Flow Cytometry method:

WBC, DIFF

RET, IPF, 2-part differential for body fluid analysis (except XN-330)

Hydrodynamic focusing DC detection method:

PLT-I (Impedance), RBC, HCT

Cyanide-free SLS-haemoglobin method:

HGB

Parameters

28 Parameters

Whole blood / Pre-dilution mode:

WBC, RBC, HGB, HCT, MCV, MCH, MCHC, PLT, RDW-SD, RDW-CV, PDW, MPV, P-LCR, PCT, NEUT#, LYMPH#, MONO#, EO#, BASO#, NEUT%, LYMPH%, MONO%, EO%, BASO%, IG#, IG%, MicroR, MacroR

Optional Parameters (Except XN-330)

14 RET & IPF Parameters

RET#, RET%, IRF, LFR, MFR, HFR, RET-He, PLT-O, IPF#, IPF, RBC-He, Delta-He, HYPO-He, HYPER-He

7 Body Fluids Parameters

Body Fluid mode

WBC-BF, RBC-BF, MN#, PMN#, MN%, PMN%, TC-BF#

Throughput

CBC :
up to 60 samples/hour

CBC + DIFF:
up to 60 samples/hour

CBC: up to 70 samples/hour
CBC + DIFF: up to 70 samples/hour
CBC + DIFF + RET: up to 35 samples/hour

Body fluid mode: up to 30 samples/hour

Sample Aspiration Volume

Whole blood mode: 25 µl
Pre-dilution mode: 70 µl

Whole blood / Low WBC mode: 25 µl
Pre-dilution mode: 70 µl
Body fluid mode: 70 µl

Data Storage

Results: 100,000 samples

Patient information: 10,000 records

QC files: 99 files / analyser

QC plots : 300 plots / file

Reagent replacement history: 5,000 records

Maintenance history: 5,000 records

Dimensions and Weight

Width : 450 mm
Depth : 460 mm
Height : 510 mm
Approx. 35 kg

Width: 450 mm
Depth: 460 mm
Height: 510 mm
Approx. 35 kg

Width: 450 mm
Depth: 460 mm
Height: 440 mm
Approx. 35 kg

Width: 450 mm
Depth: 660 mm
Height: 450 mm
Approx. 53 kg

Monitor

Width: 267 mm
Depth: 205 mm
Height: 240 mm
Approx. 3 kg

References :-

1. Sang Ook Ha, Sang Hyuk Park, So Hee Park, Jae Seok Park, Jin Won Huh, Chae-Man Lim, Younsuck Koh, Sang-Bum Hong & Seongsoo Jang (2015) Fraction of immature granulocytes reflects severity but not mortality in sepsis, *Scand J Clin Lab Invest*, 75:1, 36-43, DOI: 10.3109/00365513.2014.965736
2. Nierhaus A, Klatte S, Linssen J, Eismann NM, Wichmann D, Hedke J, et al. Revisiting the white blood cell count: immature granulocytes count as a diagnostic marker to discriminate between SIRS and sepsis - a prospective, observational study. *BMC Immunol* 2013;14(1):8.
3. Fleming C, Russcher H, Lindemans J, De Jonge R. Clinical relevance and contemporary methods for counting blood cells in body fluids suspected of inflammatory disease. Vol. 53, *Clinical Chemistry and Laboratory Medicine*. 2015. p. 1689-706.
4. Seghezzi M, Buoro S, Manenti B, Mecca T, Ferrari R, Zappalà G, et al. Optimization of Cellular analysis of Synovial Fluids by optical microscopy and automated count using the Sysmex XN Body Fluid Mode. *Clin Chim Acta*. 2016;462:41-8.
5. Weimann A, Cremer M, Hernáiz-Driever P, Zimmermann M. Delta-He, Ret-He and a new diagnostic plot for differential diagnosis and therapy monitoring of patients suffering from various disease-specific types of anemia. *Clin Lab*. 2016;62(4):667-77.
6. Dadu T, Sehgal K, Joshi M, Khodaiji S. Evaluation of the immature platelet fraction as an indicator of platelet recovery in dengue patients. *Int J Lab Hematol*. 2014;36(5):499-504.
7. Sakuragi M, Hayashi S, Maruyama M, Kabutomori O, Kiyokawa T, Nagamine K, et al. Clinical significance of IPF% or RP% measurement in distinguishing primary immune thrombocytopenia from aplastic thrombocytopenic disorders. *Int J Hematol*. 2015;101(4):369-75.
8. van der Linden N, Klinkenberg LJJ, Meex SJR, Beckers EAM, de Wit NCJ, Prinzen L. Immature platelet fraction measured on the Sysmex XN hemocytometer predicts thrombopoietic recovery after autologous stem cell transplantation. *Eur J Haematol*. 2014;93(2):150-6.



Registered under Act 737
Registration No:
IVDB17531114618
IVDB8451240217

Sysmex Asia Pacific Pte Ltd
Tel +65 6221-3629 Fax +65 6221-3687
www.sysmex-ap.com

Sysmex India Pvt. Ltd
Tel +91 (22) 6112-6666
www.sysmex.co.in

PT Sysmex Indonesia
Tel +62 (21) 3002-6688 Fax +62 (21) 3002-6699
www.sysmex.co.id

Sysmex (Malaysia) Sdn Bhd
Tel +60 (3) 5637-1788 Fax +60 (3) 5637-1688
www.sysmex.com.my

Sysmex New Zealand Ltd
Tel +64 (9) 630-3554
www.sysmex.co.nz

Sysmex (Thailand) Co., Ltd
Tel +66 (2) 032-2536 Fax +66 (2) 116-5396
www.sysmex.co.th

Sysmex Philippines Inc.
Tel +63 (2) 621-2460 Fax +63 (2) 621-2432
www.sysmex.com.ph

Sysmex Vietnam Co., Ltd
Tel +84 (8) 3997-9400 Fax +84 (8) 3997-9405
www.sysmex.com.vn

Sysmex Vietnam Co., Ltd (Hanoi Branch)
Tel +84 (4) 3776-7020 Fax +84 (4) 3776-7022
www.sysmex.com.vn

Sysmex Australia Pty Ltd
Tel +61 (2) 9016-3040
www.sysmex.com.au