

Semi-Automated Urine Chemistry Analyser

UC-1000

Every Drop Counts





Technology	Reflectance photometry with CMOS camera
Test strips	MEDITAPE UC-10S, MEDITAPE UC-12S
Throughput	max. 480 samples/hour
Parameters	<p><u>Test strip parameters</u> pH, SG, GLU, PRO, BLD, URO, KET, BIL, NIT, LEU, CRE (12S only), ALB (12S only)</p> <p><u>System parameters</u> ALB/CRE ratio (12S only), PRO/CRE ratio (12S only)</p>
Data storage	5000 samples / 5000 images of test strips
Dimension / Weight	approx. 315 x 135 x 215 / 2.6kg
Power Consumption	100 to 240 V AC (50/60 Hz) 30 VA

Key benefits

Fits and integrates anywhere

- Small footprint and lightweight
- Connects to the U-WAM as part of the UN-Series
- Equally perfect as a stand-alone solution

Fast turnaround time

- Analysis takes less than a minute
- Load another test strip every 8 seconds

Convenient use

- Simple two-step 'dip & place'
- User-friendly software for intuitive operation
- Hygienic with easy cleaning procedure

Clinical values

- Albumin and creatinine testing available
- Distinguishes between red blood cells and free haemoglobin

Registered under Act 737
Registration No: IVDB68177290818

Sysmex Asia Pacific Pte Ltd
Tel +65 6221-3629 Fax +65 6221-3687
www.sysmex-ap.com

Sysmex India Pvt. Ltd
Tel +91 (22) 6112-6666
www.sysmex.co.in

PT Sysmex Indonesia
Tel +62 (21) 3002-6688 Fax +62 (21) 3002-6699
www.sysmex.co.id

Sysmex (Malaysia) Sdn Bhd
Tel +60 (3) 5637-1788 Fax +60 (3) 5637-1688
www.sysmex.com.my

Sysmex New Zealand Ltd
Tel +64 (9) 630-3554
www.sysmex.co.nz

Sysmex (Thailand) Co., Ltd
Tel +66 (2) 032-2536 Fax +66 (2) 116-5396
www.sysmex.co.th

Sysmex Philippines Inc.
Tel +63 (2) 621-2460 Fax +63 (2) 621-2432
www.sysmex.com.ph

Sysmex Vietnam Co., Ltd
Tel +84 (028) 3997 9400 Fax +84 (028) 3997 9405
www.sysmex.com.vn

Sysmex Vietnam Co., Ltd (Hanoi Branch)
Tel +84 (024) 3776-7020 Fax +84 (024) 3776 7022
www.sysmex.com.vn

Sysmex Australia Pty Ltd
Tel +61 (2) 9016-3040
www.sysmex.com.au